

EC591 Lab 6 Report Notes

Examples of measured spectra are shown at the bottom of this document

6.1A) The spectrum of light from a fluorescent source consists of multiple lines at different wavelengths (different colors) across the visible range, which add up to produce white light.

The boxcar parameter is the number of adjacent pixels in the CCD array whose signals are averaged to produce the recorded spectrum

The larger the fiber-core diameter, the larger the number of modes guided by the fiber, and therefore the larger the range of angles of the light incident on the diffraction grating. As a result, the larger the diameter, the larger the number of pixels illuminated by each wavelength.

⇒ increasing the boxcar setting or the fiber-core diameter has the effect of broadening narrow features in the recorded spectrum

6.1B The flashlight spectrum extends over a broad range across the visible and therefore again appears as white

6.1D Regardless of fiber-core diameter, the spectral resolution of this spectrometer is not sufficient to accurately measure the linewidth of typical laser light

6.2 The reference spectrum is $I_{out}^{ref}(\lambda) = \eta I_{in}(\lambda) e^{-\alpha_w(\lambda)L}$

η : coupling efficiency from the source to the fiber through the cuvette

$I_{in}(\lambda)$: spectrum of the light from the source

$\alpha_w(\lambda)$: water absorption spectrum

$L = 1 \text{ cm}$: lateral size of the cuvette

The spectrum of the light transmitted through each salt solution is

$$I_{out}(\lambda) = \eta I_{in}(\lambda) e^{-[\alpha_w(\lambda) + \alpha_s(\lambda)]L}$$

where $\alpha_s(\lambda) = C \epsilon(\lambda)$

$\epsilon(\lambda)$ \rightarrow salt extinction coefficient $[\text{cm}^{-1}/(\text{mole/liter})]$
 C \rightarrow salt concentration $[\text{moles/liter}]$
 $\epsilon(\lambda)$ \rightarrow salt absorption coefficient $[\text{cm}^{-1}]$

The transmission spectrum recorded in each measurement is

$$T(\lambda) = \frac{I_{out}(\lambda)}{I_{ref}(\lambda)} = e^{-\alpha_s(\lambda)L}$$

$$\Rightarrow \alpha_s(\lambda) = \frac{1}{L} \ln \frac{1}{T(\lambda)} = C \epsilon(\lambda)$$

↳ NOTE: C is equal to the mass concentration
[in g/liter] divided by the formula
weight [in g/mole]

With these formulas, you can fill in all the missing data in the table on page 5 of the handout. A completed table (based on old measurement results) is included at the bottom of this document.

Mystery solution:

The absorption spectrum of the mystery solution (included below) looks like a superposition of the absorption spectra of solutions A (CuSO_4) and C (CoSO_4)

At $\lambda = 800 \text{ nm}$, the absorption coefficient of solution C is negligibly small. Therefore:

$$\alpha_{\text{mix}}(\lambda = 800 \text{ nm}) \cong x_A C_A \epsilon_A(\lambda = 800 \text{ nm}),$$

where C_A is the concentration of CuSO_4 in solution A,
 x_A is the concentration of solution A in the mystery mix,
and $\epsilon_A(\lambda)$ is the extinction spectrum of CuSO_4

By the same argument, $a_{\text{mix}}(\lambda=500\text{ nm}) \cong x_c C_c \epsilon_c(\lambda=500\text{ nm})$

From the measured values of a_{mix} at 800 & 500 nm and the relevant data in the table, you can compute x_A & x_c

\Rightarrow the mystery solution is a mixture of solutions A & C with relative concentrations x_A & x_c (nominally 40% & 60%, respectively)

6.3 A The absorption spectrum of the green dye (included below) has a minimum at about 500 nm (in the green portion of the visible range), consistent with the dye color.

6.3 B Here $I_{\text{out}}^{\text{ref}}(\lambda) = \eta_1 I_{\text{in}}(\lambda)$

\hookrightarrow coupling efficiency from the source to the fiber after reflection from the white paper

$I_{\text{out}}^{\text{milk}}(\lambda) = \eta_2 I_{\text{in}}(\lambda) e^{-\mu_{\text{milk}}(\lambda)L}$ $\left. \begin{array}{l} \text{neglecting the water absorption} \\ \text{which is small} \end{array} \right\}$

\hookrightarrow spectrum of the light transmitted through the milk cuvette

η_2 : coupling efficiency from the source to the fiber through the cuvette on the plastic box

$\mu_{\text{milk}}(\lambda)$: scattering coefficient spectrum of milk

The transmission spectrum recorded in this measurement is

$$T_{\text{milk}}(\lambda) = \frac{I_{\text{out}}^{\text{milk}}(\lambda)}{I_{\text{out}}^{\text{ref}}(\lambda)} = \frac{\eta_2}{\eta_1} e^{-\mu_{\text{milk}}(\lambda)L}$$

$$\Rightarrow \mu_{\text{milk}}(\lambda) = \frac{1}{L} \ln \frac{1}{T(\lambda)} + \underbrace{\frac{1}{L} \ln \frac{\eta_2}{\eta_1}}_{\text{unknown additive constant}}$$

The resulting scattering spectrum (included below) decreases with increasing wavelength, similar to the behavior of Rayleigh scattering


6.3C Here the recorded transmission spectrum is

$$T_{\text{milk+dye}}(\lambda) = \frac{\eta_2}{\eta_1} e^{-\mu_{\text{milk}}(\lambda)L} e^{-\alpha_{\text{dye}}(\lambda)L} = T_{\text{milk}}(\lambda) e^{-\alpha_{\text{dye}}(\lambda)L}$$


$$\Rightarrow \alpha_{\text{dye}}(\lambda) = \frac{1}{L} \ln \frac{T_{\text{milk}}(\lambda)}{T_{\text{milk+dye}}(\lambda)}$$

The dye absorption coefficient obtained from this measurement (assuming $L=1\text{cm}$) will be larger than in 6.3A, because here light travels a longer average distance from the source to the fiber due to multiple scattering events


Spectrum of Fluorescent Light Source




Spectrum of Flash Light




Spectrum of LED




Spectrum of Laser Source






Absorption Spectrum of Solution A



Absorption Spectrum of Solution B





	Name	Formula	Formula weight	Mass conc. [g/liter]	Molar conc. [moles/liter]	ϵ [$\text{cm}^{-1}/(\text{m/l})$] @ $\lambda(\text{nm})$
A	Copper sulfate	$\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$	249.7	18	0.072	29.6 @ 800
B	Nickel sulfate	$\text{NiSO}_4 \cdot 6\text{H}_2\text{O}$	262.8	20	0.076	5.1 @ 700
C	Cobalt sulfate	$\text{CoSO}_4 \cdot 7\text{H}_2\text{O}$	281.1	47	0.167	11.3 @ 500
D	Chromium potassium sulfate	$\text{CrK}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$	499.4	10	0.020	35.8 @ 600
E	Potassium chromate	K_2CrO_4	194.2	0.99	0.0051	138 @ 465

