


EC591 Lab 3 Report Notes

3.1 A



mirror reflection: $M_r = \begin{bmatrix} 1 & 0 \\ 2/R & 1 \end{bmatrix}$

free-space propagation
between the two mirrors $M_p = \begin{bmatrix} 1 & d \\ 0 & 1 \end{bmatrix}$

symmetrical confocal configuration: $R = -d$
↳ because the mirrors are concave

Starting from point P in the figure above, the ray-transfer matrix

for one pass in the cavity is $M_1 = \begin{bmatrix} 1 & 0 \\ -2/d & 1 \end{bmatrix} \begin{bmatrix} 1 & d \\ 0 & 1 \end{bmatrix} = \begin{bmatrix} 1 & d \\ -2/d & -1 \end{bmatrix}$

For two passes: $M_2 = M_1^2 = \begin{bmatrix} -1 & 0 \\ 0 & -1 \end{bmatrix}$

For four passes: $M_4 = M_2^2 = \begin{bmatrix} 1 & 0 \\ 0 & 1 \end{bmatrix}$ } the distance traveled over which
a ray reproduces itself is
 $L_{tot} = 4d$

In general the free spectral range V_{FSR} is equal to c/L_{tot}

$$\Rightarrow V_{FSR} = \frac{c}{4d} \text{ in a symmetric confocal cavity}$$

3.1B

piezoelectric responsivity $a = \frac{\Delta d}{100 \times V_{pp}}$

Δd → maximum displacement of the movable mirror

V_{pp} → peak-to-peak voltage of the scan monitor signal on the oscilloscope

→ attenuation factor of the scan monitor output (as discussed in the lab handout)


Ex for $\Delta d = 480 \text{ nm}$ and $V_{pp} = 5.8 \text{ V}$: $a = 0.83 \frac{\text{nm}}{\text{V}}$

3.1C mirror spacing $d = \frac{c}{4 \times F_{\text{FSR}}}$ { $= 7.5 \text{ mm}$ for $V_{\text{FSR}} = 10 \text{ GHz}$

FWHM of transmission peaks $\delta V_{\text{FWHM}} = \frac{V_{\text{FSR}}}{F}$ { $= 50 \text{ MHz}$ for $V_{\text{FSR}} = 10 \text{ GHz}$ and $F = 200$

As the movable mirror is displaced, the cavity length changes by a negligible amount (480 nm) compared to its initial value (7.5 mm) and therefore the FSR does not change by any appreciable amount

3.2 The figure below shows a representative oscilloscope trace of the spectrum of a HeNe laser measured with a cavity with $V_{\text{FSR}} = 10 \text{ GHz}$



Two consecutive replicas of the laser spectrum (generated by two consecutive transmission peaks of the cavity, which differ in frequency by $V_{FSR} = 10 \text{ GHz}$ are separated by a time interval $\Delta T \approx 3.4 \text{ ms}$

\Rightarrow any time increment Δt in the oscilloscope trace corresponds to a

frequency increment
$$\Delta V = \frac{V_{FSR}}{\Delta T} \Delta t \approx 2.94 \frac{\text{GHz}}{\text{ms}} \times \Delta t$$

3.2A) FSR of the laser cavity $\left. \vphantom{\begin{matrix} \text{FSR of the} \\ \text{laser cavity} \end{matrix}} \right\} V_{FSR}^{\text{laser}} = \left(\frac{V_{FSR}}{\Delta T} \right) \Delta t_{\text{adjacent lines}}$

\hookrightarrow time separation between two adjacent lines in the oscilloscope trace

In the figure above, $\Delta t_{\text{adjacent lines}} \approx 0.2 \text{ ms} \Rightarrow V_{FSR}^{\text{laser}} = 590 \text{ MHz}$

If the laser cavity has flat mirrors, its length is
$$l = \frac{c}{2V_{FSR}^{\text{laser}}} \left\{ \approx 25 \text{ cm} \right.$$

3.2B) finesse of the laser cavity $\left\{ F_{\text{laser}} = \frac{V_{\text{FSR}}^{\text{laser}}}{\delta V_{\text{FWHM}}^{\text{laser}}} \right\} = 590 \text{ if } V_{\text{FSR}}^{\text{laser}} = 590 \text{ MHz}$
 and $\delta V_{\text{FWHM}}^{\text{laser}} = 1 \text{ MHz}$

3.2C) For the spectrum shown in the figure above, the FWHM of each line was measured to be

$$\delta t_{\text{FWHM}}^{\text{laser}} \approx 20 \mu\text{s} \Rightarrow \left. \delta V_{\text{FWHM}}^{\text{laser}} \right|_{\text{measured}} = \frac{V_{\text{FSR}}}{\Delta T} \delta t_{\text{FWHM}}^{\text{laser}} \approx 59 \text{ MHz}$$

This value makes sense because the measured peaks are a convolution of the laser lines with the transmission peaks of the cavity.

↓
FWHM $\delta V_{\text{FWHM}}^{\text{laser}} = 1 \text{ MHz}$

↓
FWHM $\delta V_{\text{FWHM}} = 50 \text{ MHz}$

Because $\delta V_{\text{FWHM}} \gg \delta V_{\text{FWHM}}^{\text{laser}}$, $\left. \delta V_{\text{FWHM}}^{\text{laser}} \right|_{\text{measured}} \approx \delta V_{\text{FWHM}}$

3.2D) As you rotate the polarizer, you should see two distinct sets of equally spaced lines appear and disappear

\Rightarrow { individual lines are polarized
 but
 the overall output is unpolarized

3.2 E Grasping the laser causes a small change in the temperature inside the cavity, and therefore in the refractive index of the laser gain medium. As a result, the optical length and therefore the FSR of the laser cavity also change, and you should see all the laser lines drift in the oscilloscope

3.2 F Reflecting the laser light back into the cavity can result in caothic behavior (leading to large fluctuations in the laser spectrum) through interference effects